

Monovalent Cations in Explosives

INTRODUCTION

Isocratic elution is used to separate the common monovalent cations [sodium, ammonia, monomethylamine (MMA), and potassium] in ammonium nitrate-based inorganic explosives. Divalent cations in the sample will be retained on the column. Minimum detection limits range from about 50 ppb for sodium to about 100 ppb for monomethylamine in a 50- μ L injection. The working range of concentrations extends to about 40 ppm for each cation.

The low detection limits make this a quick, direct method for the analysis of bomb residues in criminal investigations. For quality control in explosives manufacturing, a simple dilution of the sample in deionized water is the only sample preparation required.

APPLICATION AREAS

- Explosives manufacturing
- Military quality control
- Forensic investigations

CONDITIONS

Sample Loop Volume: 50 μ L
 Columns: IonPac[®] CS3 Analytical
 IonPac CG3 Guard
 Eluent: 25 mM HCl, 0.1 mM DAP•HCl
 4% Acetonitrile
 Eluent Flow Rate: 1.0 mL/min
 Suppressor: Cation MicroMembrane[™]
 Suppressor (CMMS)
 Regenerant: 100 mM Tetrabutylammonium
 Hydroxide

Regenerant Flow Rate: 10 mL/min using AutoRegen[®]
 Accessory

Detector: Conductivity, 30 μ S full scale

Injection Volume: 50 μ L




Figure 1. Common monovalent cations in inorganic explosives

PREPARATION OF SOLUTIONS AND REAGENTS

Eluent: 25 mM HCl, 0.1 DAP•HCl 4% Acetonitrile

Dissolve 0.14 g of DAP•HCl (diaminopropionic acid monohydrochloride), 25 mL of 1 M HCl stock solution and 40 mL of acetonitrile in deionized water to 1.0 L.

Regenerant: 100 mM Tetrabutylammonium Hydroxide

Dilute 200 mL of 55% TBAOH (tetrabutylammonium hydroxide) to 4 L with deionized water.

EQUIPMENT

Dionex ion chromatograph equipped with a conductivity detector.

NOTE

Small disturbances in the baseline may be due to the elution of divalent cations. Divalent cations may be cleaned from the column with an eluent consisting of 0.5 M HCl, 0.1 mM DAP•HCl. After cleaning the column with this eluent for 15 minutes, re-equilibrate the column with the recommended eluent for this application for at least 15 minutes.



MicroMembrane is a trademark and AutoRegen and IonPac are registered trademarks of Dionex Corporation.

Printed on recycled and recyclable paper.

Dionex Corporation
1228 Titan Way
P.O. Box 3603
Sunnyvale, CA
94088-3603
(408) 737-0700

Dionex Corporation
Salt Lake City Technical Center
1515 West 2200 South, Suite A
Salt Lake City, UT
84119-1484
(801) 972-9292

Dionex U.S. Regional Offices
Sunnyvale, CA (408) 737-8522
Westmont, IL (630) 789-3660
Houston, TX (281) 847-5652
Atlanta, GA (770) 432-8100
Marlton, NJ (856) 596-0600

Dionex International Subsidiaries
Austria (01) 616 51 25 Belgium (015) 203800 Canada (905) 844-9650 China (852) 2428 3282 Denmark (45) 36 36 90 90
France 01 39 30 01 10 Germany (06126) 991-0 Italy (06) 66 51 50 52 Japan (06) 6885-1213 The Netherlands (0161) 43 43 03
Switzerland (062) 205 99 66 United Kingdom (01276) 691722
* Designed, developed, and manufactured under an NSAI registered ISO 9001 Quality System.

