



PRODUCT MANUAL

IONPAC® CTC-1 CATION TRAP COLUMN (CTC-1, P/N 040192)

IONPAC® CTC CATION TRAP COLUMN (CTC 2mm, P/N 043132)

QUICKSTART STEPS AND LINKS

Click blue text below to get started.

1. See [Section 2, "Installation"](#) for installation and operating instructions.

©DIONEX Corporation, 1997, 1998, 2002

Document No. 034536

Revision 05

7 November 2002




TABLE OF CONTENTS

SECTION 1 - INTRODUCTION	3
SECTION 2 - INSTALLATION	4
2.1 Chemicals Required	4
2.2 Solutions Required	5
2.3 Installation of the Cation Trap Column, CTC	5
2.3.1 Equilibration of the CTC to DAP·HCl Eluents	5
2.3.2 Regeneration of the Cation Trap Column When Using DAP·HCl Eluents	6
2.3.3 Regeneration and Equilibration of the CTC to Hydrochloric Acid, Sulfuric Acid or Methanesulfonic Acid	6
2.3.4 Initial CTC-1 Start-up Procedure for Trapping Ammonia in Amino Acid Eluents	7
SECTION 3 - TROUBLESHOOTING GUIDE	8
3.1 High Back Pressure	9
3.2 Unstable Retention Times	9
3.3 High Background Conductivity	10

SECTION 1 - INTRODUCTION

The IonPac® Cation Trap Column (CTC-1 for 4-mm operation, P/N 040192) and the IonPac Cation Trap Column (CTC 2-mm for 2-mm operation, P/N 043132) strips trace cationic contaminants out of the eluent and prevents them from reaching the guard and analytical columns. The cationic contaminants which are removed from the eluent by the CTC often interfere with the precision of trace cation determinations. Its use is highly recommended for high sensitivity gradient cation analysis. When performing a cation exchange application that involves an eluent gradient, a CTC should be installed between the gradient pump and the injection valve.

Table 1
IonPac CTC-1/CTC (2-mm) Packing Specifications

Column	Particle Diameter µm	Substrate ^a X-linking %	Latex Diameter nm	Latex X-Linking %	Column Capacity meq/column	Functional Group	Hydrophobicity
CTC-1 9 x 24 mm	500	8	N/A	N/A	3.0	Sulfonic acid	Low
CTC (2-mm) 4 x 35 mm	500	8	N/A	N/A	0.8	Sulfonic acid	Low

^a microporous divinylbenzene/styrene polymer

Always remember that assistance is available for any problem that may be encountered during the shipment or operation of DIONEX instrumentation and columns through the DIONEX North America Technical Call Center at 1-800-DIONEX-0 (1-800-346-6390) or through any of the DIONEX Offices listed in, "DIONEX Worldwide Offices."

SECTION 2 - INSTALLATION

The IonPac Cation Trap Column, CTC is filled with high capacity cation exchange resin. The primary application of the CTC is to strip cationic contaminants, such as metallic cations, from the acid eluents used in applications requiring the gradient elution of cations. These contaminants often interfere with the precision of trace level analysis.

The CTC-1 can also be used as an eluent ammonia trap in cationic amino acid analyses (see Section 2.3.4, "CTC-1 Initial Start-up for Trapping Ammonia in Amino Acid Eluents").

The CTC is installed in place of the high pressure Gradient Mixer that is normally positioned between the gradient pump pressure transducer and the injection valve.

2.1 Chemicals Required

Make sure that all eluents are made with high purity chemicals. Reagent grade inorganic chemicals should always be used to prepare ionic eluents. Whenever possible, inorganic chemicals that meet or surpass the latest American Chemical Society standard for purity should be used. These inorganic chemicals will detail the purity by having an actual lot analysis on each label.

Table 2
Recommended Eluents for Cation Analyses

P/N	Analytical Column	Eluent Type	Cation Trap Column
030831	HPIC™ CS1	HCl ¹	CTC-1
035371	IonPac CS2	HCl	CTC-1
037024	IonPac CS3	DAP·HCl	CTC-1
037028	IonPac CS5	Oxalic acid or PDCA ⁴	CTC Not Recommended
039591	IonPac Fast Cation-I	DAP·HCl ³	CTC-1
039633	IonPac Fast Cation-II	DAP·HCl	CTC-1
043015	IonPac CS10	DAP·HCl	CTC-1
043127	IonPac CS11	DAP·HCl	CTC (2-mm)
044001	IonPac CS12 (4-mm)	HCl or MSA ² or H ₂ SO ₄ ¹	CTC-1
044019	IonPac CS12 (2-mm)	HCl or MSA or H ₂ SO ₄	CTC (2-mm)
046073	IonPac CS12A (4-mm)	HCl or MSA or H ₂ SO ₄	CTC-1
046075	IonPac CS12A (2-mm)	HCl or MSA or H ₂ SO ₄	CTC (2-mm)
044123	IonPac CS14 (4-mm)	HCl or MSA or H ₂ SO ₄	CTC-1
044121	IonPac CS14 (2-mm)	HCl or MSA or H ₂ SO ₄	CTC (2-mm)

Notes

- #1: Use only concentrated ULTREX® grade or BAKER INSTRA-ANALYZED® for trace metals H₂SO₄ or HCl.
- #2: Use only reagent grade MSA (methanesulfonic acid).
- #3: Use DIONEX DAP (DL-2,3-diaminopropionic acid, DAP·HCl, P/N 039670) Reagent for cation eluents.
- #4: Oxalic acid and PDCA (Pyridine-2,6-dicarboxylic Acid, PDCA, P/N 039671) will strip metals from the CTC.

The deionized water used to prepare eluents should be **degassed Type I Reagent Grade Water** having a specific resistance of 18.2 megohm-cm. The deionized water should be free of ionized impurities, organics, microorganisms and particulate matter larger than 0.2 µm. Bottled HPLC-Grade Water should not be used since most bottled water contains an unacceptable level of ionic impurities. Finally, degas all deionized water prior to preparing any eluents.

2.2 Solutions Required

WARNING

Hydrochloric acid (HCl) vapors are very corrosive. Avoid breathing the vapors. Dilutions of HCl from the concentrated acid (38%) should be made in a fume hood.

Calculate the amount (in grams) of concentrated acid that you need to add to a 1-liter volumetric flask. See Table 3, "Acid Stock Solution Formulations." For example, if the HCl concentration is 38%, you need to weigh out 95.95 grams of concentrated HCl to obtain a 1.0 M HCl solution. Carefully add this amount of HCl to a 1-liter volumetric flask containing about 500 mL of deionized water with a specific resistance of 18.2 megohm-cm. Then dilute the solution to the 1-liter mark and mix thoroughly.

Table 3
Acid Stock Solution Formulations

Stock Solution Concentration	Type	Component		g	Final Volume L
		MW	%		
1.0N	HCl	36.47	38	95.95	1
See Note 1	DAPHCl	140.57	100	N/A	N/A
1.0N	H ₂ SO ₄	98.08	98	50.04	1
1.0N	MSA	96.10	100	96.10	1

Note 1

The CTC conversion solution is 50 mN HCl/20 mM DL-2,3-diaminopropionic acid monohydrochloride (DAP). See Section 2.3.1, "Equilibration of the CTC to DAP·HCl Eluents" for details of conversion. DAP·HCl is an expensive reagent. Assess your needs before making large quantities of stock solution.

2.3 Installation of the Cation Trap Column, CTC

- A. Remove the high pressure Gradient Mixer installed between the gradient pump pressure transducer and the injection valve.
- B. Connect the line from the gradient pump pressure transducer to the inlet of the CTC. DIONEX trap columns are packed with low efficiency resin because they are not positioned in the analytical pathway (injection valve--guard and analytical columns--suppressor--detector). It is not important which end of the trap column is initially designated as the inlet or outlet end of the column but after the trap column is installed in the Ion Chromatograph, it is wise not to reverse the column because the inlet end of the column will concentrate both particulates and ionic eluent contaminants. With this in mind, DIONEX places a flow direction arrow on the label of the CTC (see Figure 1, "IonPac Cation Trap Hardware Configurations").
- C. Connect a short length of liquid line from the outlet end of the CTC and direct it to a waste container.

2.3.1 Equilibration of the CTC to DAP·HCl Eluents

- A. Prepare a 500 mL solution of 50 mN HCl/20 mM DL-2,3-diaminopropionic acid monohydrochloride (DAP). DAP can be obtained from DIONEX (P/N 039670).
 1. Weigh out 1.4 g DL-2,3-diaminopropionic acid monohydrochloride (DAPHCl, MW 140.57) into a 500 mL volumetric flask.
 2. Add 25 mL of 1 N HCl stock solution to the flask.

3. Use degassed, deionized water having a specific resistance of 18.2 megohm-cm to dilute the eluent to 500 mL. Mix thoroughly to dissolve.
- B. In the following conversion procedure it is assumed that the CTC-1 is operated at 2 mL/min and the CTC (2-mm) is operated at 0.5 mL/min. As a general rule, the flow rates and eluent volumes used on the CTC (2-mm) are approximately 1/4 of the flow rates and eluent volumes used on the CTC-1.

The initial conversion of the CTC to the DAP·HCl form is achieved by pumping the conversion solution through the CTC-1 or CTC (2-mm) for approximately 4 hours or overnight at 1/4 of the above flow rates. In subsequent regenerations of the CTC after extended system use, the conversion solution through the CTC at the above flow rates for 30 minutes. This 30 minute regeneration is adequate to remove the contaminants retained during use.

- C. Prepare the appropriate DAP·HCl eluent for your application (consult the instructions shipped in the analytical column Product Manual). Pump the eluent described for the application through the CTC for 30 - 60 minutes. The column should be equilibrated within 30 - 60 minutes.
- D. Disconnect the waste line from the CTC and insert it into the ion chromatograph eluent flow path by connecting the outlet of the CTC to port #1 of the injection valve using a short length of tubing with an ID of 0.010" or less.
- E. Turn on the pump. Wait until the baseline has stabilized. Run two standard calibration runs. If the successive injections show retention times for a given solute within 5%, the CTC is fully converted and equilibrated. You may not see significant differences in background conductivity with the addition of the CTC to the system.

2.3.2 Regeneration of the Cation Trap Column When Using DAP·HCl Eluents

To regenerate the CTC when using DAP·HCl columns, repeat the conversion steps in Section 2.3.1, "Equilibration of the CTC to DAP·HCl Eluents," Steps A-E, but pump 60 mL rather than 450 mL of the conversion solution prepared in Step A through the CTC at 2 mL/min. Since the column is already in the DAP·HCl form, this 30 minute regeneration is adequate to remove the contaminants collected during previous analysis.

2.3.3 Regeneration and Equilibration of the CTC to Hydrochloric Acid, Sulfuric Acid or Methanesulfonic Acid

This procedure is used for the general removal of contaminants from the CTC. It is also used to remove DAP from the CTC when converting it for use with HCl, MSA or H₂SO₄ eluents.

In the following conversion procedure it is assumed that the CTC-1 is operated at 2 mL/min and the CTC (2-mm) is operated at 0.5 mL/min. As a general rule, the flow rates and eluent volumes used on the CTC (2-mm) are approximately 1/4 of the flow rates and eluent volumes used on the CTC-1.

- A. Disconnect the guard and analytical columns from the injection valve. Install the CTC after the injection valve. Collect the effluent from the CTC in a waste container.
- B. Dilute the acid stock solution to obtain a 100 mN acid concentration of the new eluent acid type (HCl, MSA or H₂SO₄).
- C. Rinse the CTC for approximately 1 hour with the 100 mN acid solution.
- D. Equilibrate the CTC for approximately 1 hour with the strongest eluent used in the application.
- E. Reconnect the guard and analytical columns to the outlet of the injection valve. Reinstall the CTC between the gradient pump pressure transducer and the injection valve (see Section 2.3, "Installation of the Cation Trap Column, CTC"). Resume operation.

2.3.4 Initial CTC-1 Start-up Procedure for Trapping Ammonia in Amino Acid Eluents

- A. Remove the cation exchange resin (in the hydrogen form) from the IonPac Cation Trap Column and boil it in 0.5 M NaOH (or 0.5 M LiOH if using a lithium buffer system) for about 1 hour. Rinse with large amounts deionized water with a specific resistance of 18 megohm-cm. Finally rinse with buffer A to equilibrate the resin.
 - B. Loosely pack the resin back into the IonPac Cation Trap Column body. Assemble a short length of 1/8 inch line. Connect one end of the line to the outlet of the IonPac Cation Trap Column and the other end to the bulkhead inlet of the low pressure side of the pump. Connect the line from the eluent reservoir to the inlet end of the IonPac Cation Trap Column. Reprime the pump.
 - C. Follow steps A and B to prepare an IonPac Cation Trap Column in the appropriate salt form for the B buffer line.
-

SECTION 3 - TROUBLESHOOTING GUIDE

The purpose of the Troubleshooting Guide is to help you solve operating problems that may arise while using IonPac CTC columns. For more information on problems that originate with the Ion Chromatograph (IC) or the suppressor, refer to the Troubleshooting Guide in the appropriate operator's manual. If you cannot solve the problem on your own, contact the nearest DIONEX Office (see, "DIONEX Worldwide Offices").

Table 4
CTC-1/CTC (2-mm) Troubleshooting Summary

Observation	Cause	Action	Reference Section
High Back Pressure	Unknown	Isolate Blocked Component	3.1
	Plugged Column Bed Supports	Replace Bed Supports	3.1 D
	Other System Components	Disconnect, Replace	Component Manual
Unstable Retention Times	Unequilibrated System	Lengthen First Eluent Time before Inject	3.2
	Bad Eluents	Remake Eluents	Section 2
High Background Conductivity	Bad Eluents	Remake Eluents	Section 2
	Contaminated Columns	Clean Column	Section 3.3
	Contaminated CSRS or CMMS	Clean Suppressor	Component Manual
	Contaminated Hardware	Clean Component	Component Manual

3.1 High Back Pressure

If the CTC is the cause of high back pressure, its inlet bed support may be contaminated. To change the bed support follow the instructions below using one of the two spare bed supports included in the Ship Kit.

- A. Disconnect the CTC from the Ion Chromatograph.

NOTE


Whenever you disconnect a liquid line with a ThermoFlare Fitting, inspect the opening of the line to be sure that it is not occluded. When reconnecting the liquid line, be careful not to over torque the bolt. Overtightening the bolt can seal off the end of the liquid line at the ThermoFlare.

- B. Using two open-end wrenches, carefully unscrew the inlet (the end with the label) column end fitting.
- C. Turn the end fitting over and tap it against a bench top or other hard, flat surface to remove the bed support assembly. If the bed support must be pried out of the end fitting, use a sharp pointed object such as a pair of tweezers, but be careful that you do not scratch the walls of the end fitting. Discard the old bed support assembly.
- D. Place a new bed support assembly into the end fitting. Carefully screw the end fitting onto the column so that the seal washer seats properly between the end fitting and the end of the column.

CAUTION

If the column tube end is not clean when it is inserted into the end fitting, particulate matter may prevent a proper seal between the end of the column tube and the bed support assembly. If this is the case, additional tightening may not seal the column but instead damage the column tube or the end fitting.

- E. Screw the end fitting onto the column until it is finger tight and then using wrenches, tighten it an additional 1/4 turn (25 in x lb). Tighten further only if leaks are observed.
- F. Reconnect the CTC to the Ion Chromatograph and resume operation.



3.2 Unstab

Figure 1
IonPac Cation Trap Hardware Configurations

Column Line Description	Assembly P/N	Component Descriptions
CTC-1	045287 039037 042772	End Fittings (1) Bed Support Assembly (2) Plug (not shown)
CTC (2-mm)	052809 042955 042772	End Fittings (1) Bed Support Assembly (2) Plug (not shown)

-
- A. If the column is not fully equilibrated after conversion, retention times for a given solute may vary more than 5% between injections. Equilibrate for another 30 minutes and retest.
 - B. If the conversion of the column is not complete retention times for a given solute may not be with 5% of one another. Repeat the conversion steps in Section 2.3.2, "Equilibration of the CTC to DAP·HCL Eluents," steps A - E.

3.3 High Background Conductivity

- A. If the CTC has become expended after extended use, calcium and/or magnesium peak artifacts will be observed during the analysis of blanks (such as deionized water). This happens when eluent contaminants normally not trapped by the CTC are concentrated on the guard or analytical column during equilibration with weak eluent. After a stepwise or gradient change to a stronger second eluent, they elute causing quantification interferences. These eluent contaminant peaks in the blank analysis may also be observed when the CTC is removed from the system. Were the eluents formulated correctly and did the chemicals used to make them have the required purity (see Section 2.1, "Chemicals Required")? It is important that hydrochloric acid eluents are always be prepared from trace metal grade concentrated hydrochloric acid.
 - B. Is the CTC installed in front of the injection valve in the analytical flow path (see Section 2.3, "Installation of the Cation Trap Column, CTC")? If the background conductivity is high when the CTC is not installed in the system between the gradient pump pressure transducer and the injection valve (in place of the high pressure gradient mixer), the eluent contains measurable cationic contaminants. If a new CTC or a freshly regenerated CTC (see Section 2.3, "Installation of the Cation Trap Column, CTC") is installed and the background conductivity decreases, the column is trapping eluent contaminants.
 - C. Temporarily replace the CTC, guard and analytical columns with two liquid lines between the pump and injection valve and injection valve and the Cation MicroMembrane Suppressor. Observe the background conductivity. Replace the CTC and observe the background conductivity. If the eluent is freshly prepared with high purity chemicals and if the background conductivity decreases when the CTC is removed, then the CTC is the source of the background conductivity and needs to be regenerated. To regenerate a CTC that is overloaded with contaminants, follow the steps in Section 2.3.3, "Regeneration of the Cation Trap Column After Using DAP·HCl Eluents").
-